

New biomarker: Ventana FOLR1 in ovarian cancer diagnostics

Advancing Tubo-Ovarian Cancer Management Beyond Chemotherapy

Emilia Andersson
MD, PhD, Pathologist
Oncology lead Roche Diagnostics Medical and Government Affairs EMEA-LATAM

Disclaimer


- Emilia Andersson is a full time employee of Roche Diagnostics
- The VENTANA FOLR1 (FOLR1-2.1) RxDx Assay is approved in the EU and Roche Diagnostics is providing this medical and scientific education to prepare pathologists to be confident in the interpretation of the VENTANA FOLR1 Assay to ensure there is no delay in identifying patients eligible for mirvetuximab soravtansine

VENTANA FOLR1 (FOLR1-2.1) RxDx Assay is a qualitative immunohistochemical assay using mouse monoclonal anti-FOLR1 clone FOLR1-2.1 intended for use in the assessment of folate receptor alpha (FOLR1) in formalin-fixed, paraffin-embedded epithelial ovarian cancer, including primary peritoneal cancer and primary Fallopian tube cancer (EOC) tissue specimens by light microscopy. The OptiView DAB IHC Detection Kit is used for staining on a BenchMark ULTRA instrument. FOLR1 expression is determined by the percentage of tumor cells (% TC) with membrane staining

at 2+ and 3+ intensity levels. VENTANA FOLR1 (FOLR1-2.1) RxDx Assay is indicated as an aid in identifying patients with EOC who may be eligible for treatment with the therapy listed in Table 1 for the indication and FOLR1 status in accordance with the approved therapeutic product labeling.

■ **This presentation is for educational purposes only**


S18 - 4341 A1 - 2



VENTANA FOLR1 (FOLR1-2.1) RxDx Assay

Disease context

Tubo-Ovarian cancer (OC) background



¹ Globocan database, <https://gco.iarc.fr/>, accessed June 2023, SEER database, <https://seer.cancer.gov/>, accessed June 2023

² SEER database: <https://seer.cancer.gov/>, accessed in Jan. 2022


³ De Leo, A, et al. What Is New on Ovarian Carcinoma: Integrated Morphologic and Molecular Analysis Following the New 2020 World Health Organization Classification of Female Genital Tumors. *Diagnostics*. 2021; 11(4):697. <https://doi.org/10.3390/diagnostics11040697>

⁴ NCCN Guidelines for Epithelial Ovarian Cancer Version 1.2022 (www.NCCN.org)

⁵ Ledermann, J. A., et al. (2013). "Newly diagnosed and relapsed epithelial ovarian carcinoma: ESMO Clinical Practice Guidelines for diagnosis, treatment and follow-up." *Ann Oncol* 24 Suppl 6: vi24-32. DOI: 10.1093/annonc/mdl333

⁶ Lheureux, S., et al. (2019). "Epithelial ovarian cancer." *Lancet* 393(10177): 1240-1253. DOI: 10.1016/S0140-6736(18)32552-2

Advanced ovarian cancer: a “chronic” disease with multiple relapses



PFI: platinum-free interval or duration of disease control without chemotherapy.


- Length of response shortens with each line of therapy
- Shorter platinum-free interval, lower chance of response, and shorter PFS
- Around 20% relapse within 6 months of 1st line therapy

Platinum is not the best option when:

- Progression during platinum
- Early symptomatic progression
- Platinum intolerance

González-Martín, A., et al. *Annals of Oncology* 34.10 (2023): 833-848.

Advancements in Treatment for Ovarian Cancer¹⁻⁴



1. Ovarian cancer progress timeline. American Society of Clinical Oncology. Accessed July 21, 2022. <https://www.asco.org/research-guidelines/cancer-progress-timeline/ovarian-cancer>

2. Drugs@FDA: FDA-approved drugs. US Food and Drug Administration. Accessed March 31, 2022. <https://www.accessdata.fda.gov/scripts/cder/daf/> 3. Lheureux S, et al. *CA Cancer J Clin.* 2019;69(4):480-304. 4. Manzano A, Ocaña A. *Cancers (Basel).* 2020;12(8):2223.


Folate receptor-alpha (FR α)

Background and therapeutic targeting in cancer

- Folate is a water soluble vitamin essential for **DNA and RNA synthesis** in normal cell growth and development

- Expression of membrane bound cell surface **FR α** is high during periods of rapid cell growth and division, such as **pregnancy and embryogenesis**, but is **limited** in most (not all) **adult tissues**, as they utilize an alternative transporter for folate uptake

- Due to FR α 's role in cell growth, FR α is **overexpressed in certain cancers**
- The combination of restricted expression in normal tissue and overexpression in tumor tissue makes FR α an ideal clinical target



Tumors shown to overexpress* FR α :
Breast, Lung, Endometrial, Ovarian



*derived from the literature using different detection techniques, numbers of specimens and thresholds for determining expression


Review of Folate receptor-alpha (FR α) targeting clinical trials

Mirvetuximab soravtansine

Mirvetuximab soravtansine-gynx (ELAHERE®)

Drug mechanism

- ELAHERE is an antibody-drug conjugate (ADC) targeting the FR α protein¹ on the cell surface and is internalized by the cell
- As the ADC is broken down, the cytotoxic drug is released, resulting in cell death including bystander effects due to diffusion of toxic metabolites
- By selectively targeting tumors express FR α , physicians can target the patient's cancer and reduce unwanted drug side effects
- **Efficacy correlates directly with FR α protein expression levels**
Identifying the level of FR α present in a patient's tumor is critical to ensure efficacy of ELAHERE



Developed by ImmunoGen, Inc. (now AbbVie)

Phase 2 & 3 MIRV Clinical Development Timeline

Overview of key clinical trials

FORWARD 1 Double arm

Did not meet primary endpoint but promising activity in a predefined subset of patients with high FOLR1 expression

SORAYA Single arm

Robust evidence for the efficacy of mirvetuximab soravtansine vs historical data

Led to accelerated conditional FDA approval Nov 2022

MIRASOL Double arm

Confirmatory study.

Led to full FDA approval Mar 2024 and EMA approval Nov 2024


All trials used the VENTANA FOLR1 (FOLR1-2.1) RxDx Assay to determine patient FRa expression

VENTANA® FOLR1 (FOLR1-2.1) RxDx Assay


The cutoff for eligibility

- FORWARD 1 trial: Negative results
 - Patient selection: ($\geq 50\%$ with membrane staining visible at 10x)

ITT Population




- MIRASOL trial: Positive results
 - Patient selection: ($\geq 75\%$ of cells with $\geq 2+$ membrane staining intensity)




The 75% cutoff for eligibility emphasize the importance of using a clinically validated assay

Including patients with low FOLR1 expression as per CDx, does not translate in benefit from treatment.



Commercially available RTU 1




False Positive

VENTANA FOLR1 Rx Dx Assay



FOLR1 Negative

Commercially available RTU 2



False Positive


VENTANA FOLR1 Rx Dx Assay



FOLR1 Negative

	False positive	False negative
Commercially available RTU 1	40%	0%
Commercially available RTU 2	26%	14%
Two analytically validated assays compared to the FOLR1 Rx Dx Assay		

VENTANA FOLR1 (FOLR1-2.1) RxDx Assay




Analytically and **clinically** validated IHC assay to detect FRα in epithelial ovarian, primary peritoneal, and primary fallopian tube cancer.


VENTANA FOLR1 (FOLR1-2.1) RxDx Assay is the first and only FDA- and CE-IVD approved IHC companion diagnostic for determining FRα protein expression in high grade EOC patients who may benefit from ELAHERE therapy.

VENTANA FOLR1 (FOLR1-2.1) RxDx Assay

Assay configuration



VENTANA FOLR1
(FOLR1-2.1) RxDx
Assay



OptiView DAB IHC
Detection Kit



Negative Control
(monoclonal)



BenchMark IHC/ISH
instruments



Interpretation Guide /
Scoring Algorithm



Mirvetuximab soravtansine

VENTANA FOLR1 (FOLR1-2.1) RxDx Assay

Specimen Types



Validated specimens


- Primary ovarian cancer, fallopian tube cancer and primary peritoneal cancer
- Formalin-fixed and paraffin embedded
- Archival or fresh tissues from resections, excisions, and biopsies
- Primary and metastatic sites

Non-validated specimens

- Cytology samples
- Decalcified metastatic bone lesions

System level control

- Fallopian tube with 2+ staining intensity in the epithelial cells



Estrogen-Mediated Repression of FOLR1

Cell culture models and clinical studies confirm that **Estrogen Receptor (ER) signaling actively represses** the FOLR1 gene.

The Biological Mechanism:

High Serum Estrogen → ER Binding to FOLR1 Promoter → Transcriptional Repression → Low Basolateral Protein Expression.

Pregnancy Effect:


2nd/3rd trimester serum estrogen is >500-fold higher than postmenopausal levels, causing profound FOLR1 repression in normal fallopian tissue.

Tamoxifen Effect:

Patients on Tamoxifen (ER antagonist) show significantly **higher** FOLR1 expression, as the estrogenic repression is pharmacologically blocked

Impact of Menopause on Basolateral IHC

Median Basolateral H-Score by Menopausal Group (n=51)



VENTANA FOLR1 (FOLR1-2.1) RxDx Assay

Specimen Pre-analytical Factors



Specimens should be:


- Fixation (12-72 hours) with 10% neutral buffered formalin
- 6 hour minimum fixation for small biopsies
- Sectioned at a thickness of 4µm and placed onto positively charged slides
- Cut slide stability: no more than 45 days

***Fixation with alcohol based fixatives is not recommended**

Fixation Time (hrs)	10% NBF
1*	
6*	
12	
24	
48	
72	

VENTANA FOLR1 (FOLR1-2.1) RxDx Assay

Staining Interpretation







VENTANA FOLR1 (FOLR1-2.1) RxDx Assay

Scoring Criteria

- **membrane staining is included in scoring**
 - **Partial and/or complete**
 - **Apical, circumferential distribution**
 - **When variable intensities are seen in a single cell, the **highest intensity** is scored**
- **Excluded from scoring:**
 - **Necrotic, crushed or cauterized cells**
 - **Cytoplasmic staining**

VENTANA® FOLR1 (FOLR1-2.1) RxDx Assay


Tumor cell intensity characteristics

0	1+	2+	3+
			
Absence of any detectable signal. Negative cases may still exhibit pale grey cytoplasmic and/or membranous discoloration.	Signal intensity is characterized by a faint gold/light brown hue	Signal intensity is characterized by a chocolate brown to thickened dark brown, black hue	


Note: Relative ease of visualization of membrane staining can be impacted by presence/absence of concurrent cytoplasmic staining and its intensity

VENTANA FOLR1 (FOLR1-2.1) RxDx Assay


Staining Characteristics





Membrane Staining



Apical Staining



Dot Like Staining



Cytoplasmic

Though rare, EOC cases can exhibit dot-like or apical staining that can make interpretation challenging. Both are considered as membranous staining.

FRa cytoplasmic staining should not be scored


VENTANA FOLR1 (FOLR1-2.1) RxDx Assay

Heterogeneity

Heterogeneous Staining:

Pay attention to non staining areas.

Risk of overscoring



VENTANA FOLR1 (FOLR1-2.1) RxDx Assay


Apical staining

Apical Staining:

When there is apical staining in a cell, the entire cell should be **considered positive**.

When variable intensities are seen in a single cell, **the highest intensity is scored**

Risk of underscoring



VENTANA FOLR1 (FOLR1-2.1) RxDx Assay

Dot-like staining pattern


Dot-like Pattern Staining:

Apical staining in gland forming morphology.

The entire cell should be **considered positive**.

When variable intensities are seen in a single cell, **the highest intensity is scored**


Risk of underscoring





VENTANA FOLR1 (FOLR1-2.1) RxDx Assay

Cytoplasmic and membrane staining

Cytoplasmic Staining:
Only include cells with
membrane staining.
Risk of overscoring



Dot-like and apical staining pattern



All cells touching the “dot” or the apical surface are included in the scoring

Positivity rate seems to be higher in primary vs metastatic tissue


New data from an Asian cohort (Chen et al. 2026) validates previous Western findings (Previs et al.), demonstrating global consistency in biomarker prevalence.

Study	Region	N	Positivity Rate*	Key Observation
Previs et al.	USA	432	36.3%	Primary tumors (44.4%) > Metastatic (32.5%)
Chen et al. (2026)	China	313	40.9%	Primary tumors (44.2%) > Metastatic (32.2%)

*Positivity defined as ≥75% of tumor cells staining at ≥2+ intensity.

1. Previs RA, et al. Analysis of real world FRa testing in ovarian, fallopian tube, and primary peritoneal cancers. *Gynecol Oncol.* 2025;192:102-110.
 2. Chen X, et al. Performance of the VENTANA FOLR1 Assay for folate receptor alpha: Real-world evidence from 313 Chinese participants. *Pathol Res Pract.* 2026;279:156359.

No significant Impact of Chemotherapy on FR α Expression



Key Finding: Treatment status does not significantly alter FR α expression. The positivity rate in post-chemotherapy samples is statistically comparable to chemo-naïve samples.

Optimizing Specimen selection

Resection vs. Biopsy


Due to tumor heterogeneity, the type of specimen collected significantly impacts the likelihood of a positive result.

- **Excision/Resection:** Highest yield.
- **Biopsy:** Significantly lower yield.
- **Fluid:** Lowest yield.

Recommendation: Prioritize excision or resection samples from primary tumor whenever possible to minimize the risk of false negatives due to sampling error.

Figure 4. FR α Prevalence Among Ovarian Cancer Samples by Acquisition Procedure^a (n=1574)

Sample acquisition procedure significantly impacted FR α prevalence, with resection showing the highest prevalence of FR α -positive expression ($P<0.001$)



^aA minimum of 100 viable neoplastic cells is recommended for FOLR1 testing.

For more cases and background

[Education.ventana.com](https://education.ventana.com)

or reach out to your Roche contact.

Thank you

for your attention

VENTANA, BENCHMARK are trademarks of Roche.

© 2024 Roche

All other trademarks are the property of their respective owners.

Roche Diagnostics Poland